

# Age Dating Groundwater



**Brad Esser**

**Lawrence Livermore National Laboratory**

**August 13, 2013**

# Age-Dating Groundwater: Talk Outline

- **Introduction**

- What is groundwater age?
- Why would anyone want to know groundwater age?
- How do you determine groundwater age?

- **INTRODUCED NOBLE GAS TRACERS:**

**Groundwater age dating for young water (0-2 years) using an added tracer.**


- *Application:* Recycled wastewater reuse; Managed aquifer recharge
- *Opportunity:* Commercialize LLNL NG-MIMS

- **TRITIUM/HELIUM-3 AGE DATING:**

**Groundwater age dating for older water (2-50 years) using intrinsic tracers**

- *Application:* Contaminant attribution; Groundwater monitoring;  
Basin characterization
- *Opportunity:* Make method more available to California water community


# What is groundwater age?



**Groundwater age is the travel time between recharge and discharge.**

Groundwater age usually refers to the residence time of water in the subsurface between its isolation from the atmosphere just below the water table to its discharge to a well, spring, river, lake or ocean.



# Why would you want to determine groundwater age?



Groundwater age dating can provide valuable constraints on **when and where groundwater was recharged**, and on the **age and source of groundwater contaminants**.

# Why would you want to determine groundwater age?

- **Vulnerability assessment**
  - Rapid survey of vulnerability: young = vulnerable  
(When will MTBE or nitrate show up in my well?)
- **Basin characterization**
  - Sustainable use or groundwater mining?
  - What are suitable areas for managed aquifer recharge
- **Contaminant attribution**
  - Is the contamination from legacy or current source?  
(Old irrigated agriculture or new residential septic?)
  - Is the contamination from local or non-local sources?
- **Groundwater monitoring**
  - Regional monitoring of dairies and irrigated lands
  - Assessment of best management practices
- **Regulatory compliance:**  
**managed aquifer recharge of recycled water**
  - Retention time regulations for drinking water wells.



# How do you determine groundwater age?




*Recharge of treated wastewater requires a subsurface residence time of six months to achieve a 6-log reduction in pathogens.*

| Method                   | General accuracy | General level of effort      | Retention time (months) | Safety factor |
|--------------------------|------------------|------------------------------|-------------------------|---------------|
| <b>Formula (Darcy's)</b> | Poor             | Some info on aquifer         | 24                      | 4.0           |
| <b>3-D model</b>         | Fair             | Lot of info on aquifer       | 12                      | 2.0           |
| <b>Intrinsic tracer</b>  | Better           | Sampling of existing tracers | 9                       | 1.5           |
| <b>Added tracer</b>      | Desired          | Track added tracer           | 6                       | 1.0           |

California Department of Public Health, 2008. Draft Groundwater Recharge Reuse Regulations (Title 22, California Code of Regulations; Division 4,. Environmental Health; Chapter 3. Recycling Criteria).

**Managed aquifer recharge of recycled wastewater requires tracing groundwater on a 0-2 year time scale.**

# Age tracers for short travel times (0-2 years): Managed aquifer recharge



1 month = 1 log removal

The use of treated wastewater used for groundwater recharge requires demonstrating a subsurface residence time of months to years.

## CDHS guidelines for 6-log removal credit

- **Intrinsic Tracer:** *Sample existing tracers.*  
Demonstrate a retention time of **9 months**
- **Introduced Tracer:** *Track an added tracer.*  
Demonstrate a retention time of **6 months**




## Current tracers

- **Intrinsic: Tritium/helium-3**  
Cannot determine ages <1-2 years
- **Introduced: Sulfur hexafluoride (SF<sub>6</sub>):**  
Greenhouse gas; being phased out by CARB
- **Introduced: Isotopic noble gases:**  
Limited and expensive analysis,  
Specialized sampling

# The old approach: Introduced isotopically-enriched noble gas tracers analyzed at LLNL by Noble Gas Mass Spectrometry in a dedicated facility

## Noble gases meet the technical specs for managed aquifer recharge tracers

- They are non-reactive in aqueous systems.
- They are non-toxic and approved for potable systems.
- They can be measured in small samples (10 mL).
- Multiple tracers can be introduced at different times or to places and be measured in a single analysis.
- They are widely available and are inexpensive enough (~1\$/1000 m<sup>3</sup>) for large-scale experiments.







# The old approach: Introduced isotopically-enriched noble gas tracers analyzed at LLNL by Noble Gas Mass Spectrometry in a dedicated facility

## Limitations of the old method

- Inefficient tracer introduction:  
*Significant loss of gas from bubbling*
- Specialized sampling protocols: *Copper tubes*
- Low-throughput: *8 samples/day*
- Near-unique facility:  
*LLNL, ~\$1M setup and highly trained staff*
- No possibility for field deployment
- Low possibility for commercialization




# The new approach: Introduced noble gas tracers analyzed on a small, simple Noble Gas Membrane Inlet Mass Spectrometer (NG-MIMS)




## The new method addresses these limitations

- Efficient tracer introduction:  
*Diffusion tubing for <5% loss*
- Standard sampling protocols: *40-mL VOA vial*
- High-through: *8-10 samples/hour*
- Benchtop instrument:  
*<\$50K in commercial off-the-shelf components*
- Strong possibility for field deployment
- Strong possibility for commercialization





# The new approach: Introduced noble gas tracers analyzed on a small, simple Noble Gas Membrane Inlet Mass Spectrometer (NG-MIMS)



Fast analyses on a new **NG-MIMS** bench-top instrument made with commercial parts

## THE NEW APPROACH HAS BEEN FIELD TESTED.

- **The new approach is inexpensive:**  
Xenon tracer studies at a similar cost to SF<sub>6</sub>,
- **Tracer introduction and sample collection are simple:**  
Un-attended introduction at near 100% efficiency  
Sample collection using standard methods
- **The new approach is sensitive:**  
A small fraction (0.1%) of surface recharge water can be detected in drinking water production wells





# NG-MIMS (Introduced xenon tracer analysis): Opportunities

- **Add groundwater xenon tracing to your bag of tools**
  - LLNL does all sample processing, analysis and data reduction
  - LLNL provides guidance on dosing and sampling
  - LLNL provides or assists with data interpretation
- **Contract with LLNL directly or indirectly**
  - Work with an analytical service provider that contracts with LLNL:
    - *Standard contract quickly put in place*
    - *LLNL has a non-exclusive contract with Hydrotrace, LLC*
  - Contract directly with LLNL:
    - *Lengthier process, non-standard contract*
    - *Act as an analytical services provider of this capability to the California water resource community*

# NG-MIMS (Introduced xenon tracer analysis): Opportunities

- **Contract with LLNL to build an instrument or to lease an instrument for an extended tracer campaign**
  - We are currently developing a contract with a non-California utility with this concept
- **Use the LLNL design to build and sell a commercial instrument**
  - LLNL has an ROI but not a patent
  - There is a lot of room for improvement
    - Place instrument in a standard cabinet
    - Ruggedize for field deployment
    - Develop software
      - Integrate instrument control and data acquisition
      - Develop a data reduction module
      - Make more use friendly





# NG-MIMS (Introduced xenon tracer analysis): Intellectual Property

- **LLNL does have a record of invention in NG-MIMS**
  - M. Singleton and A. Visser, 2010, Noble Gas Membrane Mass Spectrometer, IL-12385
  - LLNL has elected not to pursue a patent at this time
- **LLNL does have reports and an accepted peer-reviewed publication**
  - **Method report** : Visser, A., Singleton, M., Hillegonds, D., Velsko, C., Moran, J.E., Esser, B.K., 2012. California GAMA Special Study: A Noble Gas Membrane Inlet Mass Spectrometry (NG-MIMS) system for water and gas samples. Lawrence Livermore National Laboratory LLNL-TR-548931, p. 22.
  - **Field study report**: Visser, A., Singleton, M., Hillegonds, D., Velsko, C., Moran, J.E., Esser, B.K., 2013. California GAMA Special Study: Rapid, Low-Cost Noble Gas Tracer Monitoring to Determine Travel Times at Recharge Operations. Lawrence Livermore National Laboratory LLNL-TR-534291, p. 40.
  - **In press manuscript**: Visser, A., Singleton, M.J., Hillegonds, D.J., Velsko, C.A., Moran, J.E., Esser, B.K., 2013. A Membrane Inlet Mass Spectrometer for Noble Gases at Natural Abundances in Gas and Water Samples, *in press*, Rapid Communications in Mass Spectrometry.
- **NG-MIMS for groundwater studies is an active field in Europe and to a lesser extent in this country**
  - **Similar but not identical instrument from Europe**: Mächler, L., Brennwald, M.S., Kipfer, R., 2012. Membrane Inlet Mass Spectrometer for the Quasi-Continuous On-Site Analysis of Dissolved Gases in Groundwater. Environmental Science & Technology 46, 8288-8296.
  - **A MIMS, but not capable of measuring xenon**: Kana, T.M., Darkangelo, C., Hunt, M.D., Oldham, J.B., Bennett, G.E., Cornwell, J.C., 1994. Membrane inlet mass spectrometer for rapid high precision determination of N<sub>2</sub>, O<sub>2</sub>, and Ar in environmental water samples. Analytical Chemistry 66, 4166-4170.

## Age tracers for longer travel times (2-60 years):


Nitrate is the most important water quality issue in California groundwater and requires age dating on a decadal time scale



*Most public-drinking water wells in California produce water with measurable tritium.*


**Groundwater monitoring for contaminant source attribution, vulnerability assessment, and basin characterization requires tracing groundwater on a 0-60 year time scale**

# GAMA Special Studies is developing new tracers to trace recent recharge







# Tritium and its decay product $^3\text{He}$ are used to determine groundwater age (<50y)



Cosmic rays



Atmospheric  
weapons testing



**Groundwater  
recharge**

Tritium DL < 1 pCi/L

Rainfall today = ~10 pCi/L

Tritium MCL = 20,000 pCi/L

# Why not just use tritium alone?

Tritium activity is affected by groundwater age, groundwater mixing *and recharge source*

| Location    | Tritium (pCi/L) | GW age (years) |
|-------------|-----------------|----------------|
| Bakersfield | 11.1            | 15             |
|             | 11.3            | 34             |
| San Jose    | 12.4            | 11             |
|             | 12.5            | 42             |

## Irrigation with pumped groundwater

- **Will not affect** tritium-helium groundwater ages
- **Will affect** tritium-only model ages

The use of the tritium/helium-3 method also gives valuable information on

- **Recharge Temperature and Elevation:**  
*Pond vs river; high vs low elevation*
- **Pre-modern recharge:**  
*Mixing between water older than and younger than 50 years*
- **Very old water components:**  
*Helium-4 age dating*
- **Excess air and style of recharge:**  
*Injection vs surface spreading*

# Tritium/helium-3 groundwater age dating and noble gas geochemistry is a signature LLNL capability

- **LLNL is one of two high-throughput noble gas mass spectrometry labs in the country**
  - 500-1,000 samples per year
  - $^3\text{He}$  accumulation for  $^3\text{H}$  determination
  - Noble gas composition by isotope dilution on a residual gas analyzer
  - $^3\text{He}/^4\text{He}$  on VG5400



*One of two 12-position sample manifolds under LabView control*

*Ate Visser, lab manager*

## Noble Gas Mass Spectrometry



The minimum detection level for  $^3\text{He}$  is about 1 zeptomole (600 atoms). We routinely measure  $10^5$  atoms of  $^3\text{He}$  with accuracy of a few percent

Over the past 15 years, LLNL has made groundwater age-dating accessible to water resource managers

# NG-MIMS (Introduced xenon tracer analysis): Opportunities

- **Add groundwater xenon tracing to your bag of tools**
  - LLNL does all sample processing, analysis and data reduction
  - LLNL provides guidance on dosing and sampling
  - LLNL provides or assists with data interpretation
- **Contract with LLNL directly or indirectly**
  - Work with an analytical service provider that contracts with LLNL:
    - *Standard contract quickly put in place*
    - *LLNL has a non-exclusive contract with Hydrotrace, LLC*
  - Contract directly with LLNL:
    - *Non-standard contract, lengthier process*
    - *Act as an analytical services provider of this capability to the California water resource community*

# NGMS (Tritium/Helium-3 Age Dating): Opportunities

- **For analytical service providers**
  - **Model 1: LLNL does all analysis**
    - Service provider directly contracts with LLNL, and bundles smaller contracts
    - LLNL analyzes tritium samples and noble gas samples
  - **Model 2: LLNL does noble gas analysis**
    - LLNL analyzes noble gas samples only
    - Service provider analyzes tritium samples by enrichment and liquid scintillation counting
  - **Model 3: LLNL does mass spectrometry**
    - Service provider degasses tritium samples;
    - LLNL analyzes degassed tritium samples, and analyzes noble gas samples

**Stretch goal: Develop a small instrument that can measure helium isotopic composition with 2% precision to allow commercial labs to provide tritium/helium-3 groundwater age dates**



**Degassing rack**



# NGMS (Tritium/Helium-3 Age Dating): Intellectual Property

- **LLNL does not have ROIs, patents or intellectual property in NGMS**
  - for determination of tritium by helium-3 accumulation, or
  - for determination of helium isotopic composition or noble gas abundance
- **LLNL offers a unique capability: High-throughput analysis of groundwater for tritium/helium groundwater age dating**
  - Past experience has shown that this is difficult to setup and maintain
    - Two high-throughput laboratories in US: LLNL and Columbia University
    - A few low-throughput university laboratories, e.g. University of Utah
  - LLNL must not compete with private industry
- **LLNL does have standard operating protocols for these analyses**
  - Visser, A., Hillegonds, D., Esser, B.K., 2013. Collection and Analysis of Groundwater for Determination of Tritium by Helium-3 Accumulation (SOP-NGMS-121 revision 5). Lawrence Livermore National Laboratory NGMS Standard Operating Procedure (LLNL-TM-623415), p. 9.
  - Visser, A., Hillegonds, D., Esser, B.K., 2013. Collection and Analysis of Groundwater for Determination of Noble Gas Abundance and Helium Isotopic Composition (SOP-NGMS-122 revision 4). Lawrence Livermore National Laboratory NGMS Standard Operating Procedure (LLNL-TM-623335), p. 15.

# Age-Dating Groundwater: Opportunities to Work with LLNL

- **INTRODUCED NOBLE GAS TRACERS:**

**Groundwater age dating for young water (0-2 years) using an added tracer.**

- *Application:* Recycled wastewater reuse; Managed aquifer recharge
- *Opportunity:* Commercialize LLNL NG-MIMS

- **TRITIUM/HELIUM-3 AGE DATING:**

**Groundwater age dating for older water (2-50 years) using intrinsic tracers**

- *Application:* Contaminant attribution; Groundwater monitoring; Basin characterization
- *Opportunity:* Make method more available to California water community

- **IF YOU ARE INTERESTED**

- **Brad Esser:** 925-422,-5247, [bkesser@llnl.gov](mailto:bkesser@llnl.gov)  
[https://www-pls.llnl.gov/?url=about\\_pls-scientific\\_staff-esser\\_b](https://www-pls.llnl.gov/?url=about_pls-scientific_staff-esser_b)
- **Aaron Tremaine:** 925-422-1284, [tremaine1@llnl.gov](mailto:tremaine1@llnl.gov)
- **State Water Board Reports Library: GAMA Special Studies**  
[http://www.waterboards.ca.gov/water\\_issues/programs/gama/report\\_depot.shtml](http://www.waterboards.ca.gov/water_issues/programs/gama/report_depot.shtml)

# Time for questions and a break

